




COVID-19 IgG/IgM Rapid Test

- Whole blood (Finger prick available), serum and plasma sample
- Easy to use - no analyser required
- Store at room temperature or refrigerator



2019-nCoV



The novel 2019 Coronavirus Pneumonia (NCP) or “COVID-19” was confirmed and named by the World Health Organisation (WHO) in January 2020 when it was discovered in Wuhan, China.

On January 30th 2020 WHO released the new Coronavirus infection pneumonia epidemic as a public emergency of international concern.

The pneumonia is caused by the novel coronavirus (SARS-CoV-2) mainly transmitted through respiratory droplets. Typical symptoms of coronavirus include fever and a cough that may progress to a severe pneumonia causing shortness of breath and breathing difficulties. Coronavirus can cause more severe symptoms in people with weakened immune systems, older people, and those with long-term conditions like diabetes, cancer and chronic lung disease. Prestige Diagnostics UK are now offering a COVID-19 IgG/IgM Rapid Test offering high accuracy and quick results.

The COVID-19 Rapid Test Cassette has been compared with a commercial PCR test, the results indicating high specificity and sensitivity: IgG results: Relative Sensitivity: 100% (CI*: 86.0% - 100%) Relative Specificity: 98.0% (95% CI*: 89.4% - 99.9%) Accuracy: 98.6% (95%CI*: 92.3% - 99.96%)


IgM results: Relative Sensitivity: 85.0% (95% CI*: 62.1% - 96.8%) Relative Specificity: 96.0% (95% CI*: 86.3% - 99.5%) Accuracy: 92.9% (95% CI*: 84.1% - 97.6%)

*Confidence Intervals



Cat No.	Production Description	Pack Size
RADCOV1	10 Min RAPID Test Device	20 Tests
RADCOV2	10 Min RAPID Test Device with Lancet and Alcohol Swab	20 Tests

VISUAL RESULTS AVAILABLE IN ONLY 10 MINUTES



AMS U.K. Limited

Ballymena Business Centre, Galgorm,
County Antrim, United Kingdom. BT42 1FL

www.amslabs.co.uk | info@amslabs.co.uk

T: +44 (0) 28 2565 6268

